


<b>Product Description</b>	<p>The Formaldehyde Loadcheck device is designed as a penetration test for use in low temperature steam / formaldehyde (LTSF) sterilizers. It is a challenge test for determining whether steam/ formaldehyde vapour is penetrating into the lumen type device.</p> <p>The Formaldehyde Loadcheck consists of a helix device in which a chemical indicator is located in the capsule at the end of a length of tube. If steam / formaldehyde vapour penetrates through to the end of the tube the indicator will show a defined colour change from blue to green.</p> <p>The Loadcheck device has a limited life, so after 200 cycles the device must be replaced.</p>											
<b>Specification</b>	<p>Colour changes from blue to completely green if steam / formaldehyde sterilant vapour penetration of the sterilant has been achieved.</p>											
<b>Applicable standards</b>	<p>No published standard for this product</p>											
<b>Product Size</b>	<p>N/A</p>											
<b>Packaging</b>	<table border="0"> <tr> <td data-bbox="389 850 519 892">Primary case contents:</td> <td data-bbox="519 850 779 892">200 indicators 1 Loadchek Device</td> <td data-bbox="779 850 876 892">Size:</td> <td data-bbox="876 850 1120 892">130mm x 110mm x 30mm</td> <td data-bbox="1120 850 1459 892">Weight: 135g Gross</td> </tr> <tr> <td data-bbox="389 934 519 976">Secondary case contents:</td> <td data-bbox="519 934 779 976">20 primary cases</td> <td data-bbox="779 934 876 976">Size:</td> <td data-bbox="876 934 1120 976">350mm x280mm x 165mm</td> <td data-bbox="1120 934 1459 976">Weight: 3kg Gross</td> </tr> </table>		Primary case contents:	200 indicators 1 Loadchek Device	Size:	130mm x 110mm x 30mm	Weight: 135g Gross	Secondary case contents:	20 primary cases	Size:	350mm x280mm x 165mm	Weight: 3kg Gross
Primary case contents:	200 indicators 1 Loadchek Device	Size:	130mm x 110mm x 30mm	Weight: 135g Gross								
Secondary case contents:	20 primary cases	Size:	350mm x280mm x 165mm	Weight: 3kg Gross								
<b>Recommended Storage</b>	<p>Store the device in cool, dry conditions before and after use and with the capsule removed to aid drying : 0 - 30°C 30 – 60 % RH</p>											
<b>Shelf Life</b>	<p>24 months from date of manufacture</p>											
<b>Active Components</b>	<p>Each indicator, which is the indicating medium, contains 0.7mg dyes and 0.7mg reagents. None of the substances used in the formulation are known carcinogens, nor do they contain any heavy metals.</p> <p>The component chemicals are potentially irritant. Being fully encapsulated, this should present no occasion for contact with the chemical components. If however, contact is made due to breakdown of the encapsulation, wash the skin with soap and water.</p> <p>If ingested wash out the mouth thoroughly and give plenty of water to drink.</p>											
<b>Safe Usage</b>	<p>The indicator and packaging contain no added heavy metals, no known carcinogens, and no added rubber latex. The indicator is considered safe when used under normal conditions.</p> <p>When used as intended, the indicator does not release any substances known to be toxic in sufficient quantities to cause a health hazard, during or after the sterilization process for which it is designated, in accordance with ISO 11140-1:2005, 5.9</p>											
<b>Disposal</b>	<p>Treat unwanted indicators and packaging in the same way as normal paper waste.</p>											